

**INDIANA DEPARTMENT OF TRANSPORTATION
MATERIALS AND TESTS DIVISION
SPENT FOUNDRY SAND TOXICITY TEST
ITM No. 215-02T**

1.0 SCOPE.

- 1.1** This test method covers the procedure for the rapid evaluation of the toxicity of spent foundry sand using luminescent marine bacterium.
- 1.2** The luminescent marine bacterium, *Vibrio fischeri*, produce light as a result of the normal biological activity for this bacterium. The test exposes the luminescent bacterium to aqueous leachate samples and measures the change in the light output at 5 and 15 minutes. Leachate samples that contain toxins cause a reduction of the light output due to a decrease in the biological activity of the bacterium. By measuring the relative change in the light output of the bacterium exposed to the leachate samples verses the change in the light output of a control sample the relative toxicity of the spent foundry sand can be determined.
- 1.3** The values stated in either acceptable English units or SI metric units are to be regarded separately as standard, as appropriate for a specification with which this ITM is used. Within the text, SI metric units are shown in parenthesis. The values stated in each system may not be exact equivalents; therefore each system shall be used independently of each other, without combining values in any way.
- 1.4** This ITM may involve hazardous materials, operations, and equipment. This ITM does not purport to address all of the safety problems associated with the ITMs use. The ITM user's responsibility is to establish appropriate safety and health practices and to determine the applicability of regulatory limitations prior to use.

2.0 REFERENCES.

2.1 ASTM Standards

- D 5660 Test Method for Assessing the Microbial Detoxification of Chemically Contaminated Water and Soil Using a Toxicity Test with Luminescent Marine Bacterium
- E 943 Terminology Relating to Biological Effects and Environmental Fate

2.2 INDIANA TEST METHODS OR PROCEDURES

- 207 Sampling Stockpiled Aggregates
- 802 Random Sampling

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom"

2.3 Other

MICROTOX ® 90% comparison test protocol, 1998 AZUR Environmental.

3.0 TERMINOLOGY. Definitions for terms and abbreviations shall be in accordance with the Department's Standard Specifications, Section 101 and ASTM E 943.

4.0 SIGNIFICANCE AND USE. The test method is used as an indicator to determine acceptance or rejection of spent foundry sand from a ferrous foundry for highway construction projects. Spent foundry sand must first conform to Indiana Administrative Code, 329 IAC 10, for Type III or Type IV restricted waste. Testing for this typing is accomplished by the foundry using accepted EPA Methods. When resampling of the spent foundry sand is required for waste reclassification, 329 IAC 10-9-4 testing per this protocol is also required.

5.0 APPARATUS. The testing apparatus shall be in accordance with ASTM D 5660.

6.0 SAMPLING. Unless otherwise directed by the Engineer, samples shall be obtained as follows:

- 6.1** Identify the volume and portion of waste foundry sand open-stockpile that will be used on construction projects. Divide this portion of the stockpile into lots of approximately 20,000 CYS (15,000 m³). Divide each lot into five sublots of equal volume.
- 6.2** Lots of less than 20,000 CYS (15,000 m³), each subplot shall be approximately 4000 CYS (3000 m³).
- 6.3** Determine random sampling locations within each subplot in accordance with ITM 802. The depth of the sampling shall never be less than 1 ft (0.3 m) or require an excavation greater than 8 ft (2.5 m).
- 6.4** Sample each subplot in accordance with ITM 207.
- 6.5** Foundry's that store spent foundry sand using a method other than an open - stockpile shall provide a QCP which provides for random sampling of each subplot of 4000 CYS (3000 m³).

7.0 SAMPLE PREPARATION. The leachate sample preparation of each subplot sample for toxicity testing shall be as follows:

- 7.1** Place 20.00 ±0.05 grams of the spent foundry sand into a flask, add 80 ±1mL of water to the flask. Cover the flask with parafilm, manually agitate to break down any large clumps, and place the flask on a shaker table at 175 to 200 RPM for 18 ±2

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX ® test kit for use as proposed in this Test Method. Use of the MICROTOX ® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX ® test kit in this manner or for any results therefrom."

h. No mechanical breaking or crushing of the sand shall be performed in order to represent actual site conditions.

7.2 After shaking allow the sample to settle, then pour the supernatant into polycarbonate centrifuge tubes and centrifuge for 16 minutes at 10000 RPM. Filter the supernatant using a 1.5 µm pore size glass fiber filter and then a 0.45 µm pore size membrane filter to remove fines. Measure the pH of the filtrate and transfer to borosilicate glass vials, cover with parafilm and cap.

7.3 Each of the leachate samples shall be tested immediately or stored at 39°F (4°C) for no more than 72 h prior to testing.

8.0 PROCEDURE. The leachate samples shall be tested in accordance with the MICROTOX ® 90% comparison test protocol in appendix A.

9.0 CALCULATION

9.1 For each leachate sample, calculate a control test mean, X_c , and a sample test mean, X_s , for each time T5 and T15.

9.2 Calculate the normalized mean difference, \bar{X} , between the means of the control test and the sample test for each time, T5 and T15, for each subplot sample.

$$\bar{X} = \frac{(X_c - X_s \times 100)}{X_c}$$

9.3 Calculate the pooled mean difference from all the samples, \bar{X}_p

$$\bar{X}_p = \frac{1}{p} \sum_{i=1}^p \bar{X}_i \quad i = 1, 2, \dots, p$$

Where p = the number of subplot samples

\bar{X}_i = the normalized mean difference of the i -th sample

9.4 If any \bar{X}_i is greater than $\bar{X} + 2\sigma$, assume $\sigma = 3.0$, then obtain a new subplot sample, include both new sample values and the original values for \bar{X}_i in the \bar{X}_p calculations.

10.0 REPORT. The report of the results shall consist of TCLP and Neutral Leachate Test results upon which the waste classification is based, certification that the ferrous spent foundry sand is Type III or Type IV as set out in the Recurring Special Provisions for use of spent foundry sand in construction projects, pH of the leachate from each sample, and the results of the MICROTOX ® 90% comparison test protocol. The MICROTOX ® results

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX ® test kit for use as proposed in this Test Method. Use of the MICROTOX ® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX ® test kit in this manner or for any results therefrom."

shall include T5 and T15 mean values for each of the sublot samples including rejected sublot samples. The dates of the sampling and testing shall be included in the test report. The test results shall be kept on file for a minimum of five years. The report shall be sent to the Environment, Planning and Engineering Division.


“The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc.”

“Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom.”

Appendix A


“The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc.”

“Strategic Diagnostics Inc. has not tested or approved the MICROTOX ® test kit for use as proposed in this Test Method. Use of the MICROTOX ® test kit as proposed herein is n no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX ® test kit n this manner or for any results therefrom.”



Page 1

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.



Page 2

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom."

Comparison Test

- Add 10.0 mL Reconstitution Solution to a test tube.
- Add 1,000 μ L OAS to the test tube, mix.
- Transfer 1,500 μ L Osmotically Adjusted Reconstitution Solution to six cuvettes. Place Osmotically Adjusted Reconstitution/cuvettes in wells:
B1, B3, B5
D2, D4
F3

Page 3

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.

Comparison Test

- Add 10.0 mL of sample to a test tube.
- Add 1,000 μ L OAS to the test tube, mix.
- Transfer 1,500 μ L Osmotically Adjusted Sample to five cuvettes. Place Osmotically Adjusted Sample/cuvettes in wells:
B2, B4
D1, D3, D5
- Wait 5 minutes.

Page 4

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom."

Comparison Test


Reagent Reconstitution

Reconstitute a vial of Microtox Acute Toxicity Reagent in the following way:

- Remove a single vial of reagent from the freezer and open it with the minimum of handling, thereby reducing warming of the vial.
- Shake and tap the vial gently to ensure the pellet of bacteria is seated on the bottom of the vial.
- Take the precooled cuvette of Reconstitution Solution from the Reagent Well, then quickly pour the solution into the opened vial.
- Swirl the vial 3 or 4 times, then quickly pour the mixture back into the cuvette and return it to the Reagent well.
- Mix the bacteria thoroughly using the pipettor by aspirating and dispensing 0.5 ml of solution at least 10 times. Reconstituted bacteria should be used within 3 hours of reconstitution. Further tests after this period require the preparation of freshly reconstituted bacteria.

Page 5

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.




Page 6

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.


"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom."



Page 7

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmni and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.



Page 8

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmni and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom."

Comparison Test

1	2	3	4	5	
Gray	Gray	Gray	Gray	Gray	A
Cyan	Magenta	Cyan	Magenta	Cyan	B
Yellow	Yellow	Yellow	Yellow	Yellow	C
Magenta	Cyan	Magenta	Cyan	Magenta	D
Yellow	Yellow	Yellow	Yellow	Yellow	E
Gray	Gray	Yellow	Gray	Gray	F

Reagent Well: Yellow

- READ zero time I_0 light levels as prompted by the computer monitor:
C1, C2, C3, C4, C5
E1, E2, E3, E4, E5
- Then immediately...

Page 9

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.

Comparison Test

1	2	3	4	5	
Gray	Gray	Gray	Gray	Gray	A
Cyan	Magenta	Cyan	Magenta	Cyan	B
Yellow	Orange	Yellow	Orange	Yellow	C
Magenta	Cyan	Magenta	Cyan	Magenta	D
Orange	Yellow	Orange	Yellow	Orange	E
Gray	Gray	Yellow	Gray	Gray	F

Reagent Well: Yellow


- ... make the following 900 μ L transfers:
B1 to C1
B2 to C2
B3 to C3
B4 to C4
B5 to C5
D1 to E1
D2 to E2
D3 to E3
D4 to E4
D5 to E5
- Touch the computer <space bar> key.

Page 10

©1998 AZUR Environmental. All rights reserved. Microtox, Mutatox and MicrotoxOmri and the AZUR logo are trademarks or registered trademarks of AZUR Environmental.

"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom."



"The trademark MICROTOX® is used with the permission of its registered owner, Strategic Diagnostics Inc."

"Strategic Diagnostics Inc. has not tested or approved the MICROTOX® test kit for use as proposed in this Test Method. Use of the MICROTOX® test kit as proposed herein is in no way sponsored, endorsed or authorized by Strategic Diagnostics Inc. and Strategic Diagnostics Inc. shall have no liability for use of the MICROTOX® test kit in this manner or for any results therefrom."